

Materials List for

# Lung Instillation and Fixation: Introducing Liquid into the Trachea and Fixing Lungs to Measure the Fixed Lung Volume

1

URL: <https://www.jove.com/video/20251>

DOI: [doi:10.3791/20251](https://doi.org/10.3791/20251)

## Materials

Name	Company	Catalog Number	Comments
Laboratory Balance	Ohaus	Adventurer Pro Model AV 313	Other balances can be used if they have a range of 1-300 g
20 g Luer intravenous catheter	Insylte		Several other possible vendors, e.g., Jelco Optiva
500 ml laboratory bottle	Various		Several other possible vendors