

Materials List for:

Isolation of Genomic DNA from Mouse Tails

Tony Zangala¹

¹Department of Physiology and Biophysics, University of California, Irvine (UCI)

Correspondence to: Tony Zangala at T.Zangala@gmail.com

URL: <https://www.jove.com/video/246>

DOI: [doi:10.3791/246](https://doi.org/10.3791/246)

Materials

Name	Type	Company	Catalog Number	Comments
STE	Buffer			100 mM Tris; pH 8.5 /5 mM EDTA /0.2% SDS /200 mM NaCl / in H2O
TE	Buffer			100 mM Tris; pH 5.0 /1 mM EDTA /10 mM NaCl