

Materials List for:

Assessment of Sensorimotor Function in Mouse Models of Parkinson's Disease

Sheila M. Fleming^{1,2}, Osunde R. Ekhator^{1,2}, Valentins Ghisays¹

¹Department of Psychology, University of Cincinnati

²Department of Neurology, University of Cincinnati

Correspondence to: Sheila M. Fleming at fleminsa@ucmail.uc.edu

URL: <https://www.jove.com/video/50303>

DOI: [doi:10.3791/50303](https://doi.org/10.3791/50303)

Materials

Name	Company	Catalog Number	Comments
Challenging beam apparatus	Starks Plastics	Segment 1= 3.5 cm Segment 2= 2.5 cm Segment 3= 1.5 cm Segment 4= 0.5 cm	Contact information: 11276 Sebring Dr. Forest Park ,OH 45240 USA Ph +1 (513) 541-4591 Fax +1 (513) 541-6773
Cylinder	Starks Plastics	15.5 cm diameter 12.7 cm height	Contact information: 11276 Sebring Dr. Forest Park ,OH 45240 USA Ph +1 (513) 541-4591 Fax +1 (513) 541-6773
Mesh Grid Wiring	Ace Hardware	1 cm ²	
Mouse Cages	Ancare	19 x 29 x 12.7 cm	
Glass	Ace Hardware	19 cm ²	
Mirror	Ace Hardware	18 x 13 cm	
Camcorder	Sony HDR-HC9		
MiniDV Tapes	Sony DVC premium	90 min long play	
Labels	AveryColor Coding Labels	6.35 mm diameter (1/4" Round)	