

Materials List for:

SpOT the Correct Tissue Every Time in Multi-tissue Blocks

Anna Coffey^{1,2}, Michael D. Johnson³, Deborah L. Berry³

¹Center for Advanced Preclinical Research (CAPR), Frederick National Laboratory for Cancer Research

²Leidos Biomedical Research, Inc.

³Department of Oncology, Georgetown University

Correspondence to: Deborah L. Berry at dlb82@georgetown.edu

URL: <https://www.jove.com/video/52868>

DOI: [doi:10.3791/52868](https://doi.org/10.3791/52868)

Materials

Name	Company	Catalog Number	Comments
Histogel Specimen Processing Gel	Thermo Scientific	HG-4000-012	http://www.thermoscientific.com/en/product/richard-allan-scientific-histogel-specimen-processing-gel.html
Tissue Marking Dye	Triangle Biomedical Sciences, Inc.	TMD-5	Any tissue marking dye would most likely be sufficient.
Arraymold Kit A 2 mm (60 core)	Arraymold	20015A	Any manual tissue arrayer would work similarly.