

Materials List for:

Assessment of Child Anthropometry in a Large Epidemiologic Study

Amy L. Louer¹, Denise N. Simon¹, Karen M. Switkowski¹, Sheryl L. Rifas-Shiman¹, Matthew W. Gillman¹, Emily Oken¹

¹Obesity Prevention Program, Department of Population Medicine, Harvard Medical School and Harvard Pilgrim Health Care Institute

Correspondence to: Emily Oken at emily_oken@harvardpilgrim.org

URL: <https://www.jove.com/video/54895>

DOI: [doi:10.3791/54895](https://doi.org/10.3791/54895)

Materials

| Name | Company | Catalog Number | Comments |
|---|-------------------------------|---|---|
| Stadiometer | Weigh and Measure, LLC | SSWM-1 | Basic Shorr board (without smooth slide features) can also be used. In order to accommodate the width of children's hips during sitting height, the base of a stadiometer should be approximately 60 cm wide or larger. |
| Bioimpedance scale | Tanita Coporation of America | TBF 300A (model is discontinued), DC-430U is comparable | Scale is used for weight and bioimpedance. Any digital, standardized scale can be used for weight only. |
| Skinfold Caliper | Holtain Limited | n/a | This model uses a dial gauge in graduations of 0.2 mm. Models with a linear gauge are also acceptable. |
| Hip/waist tape measure | Gulick II Plus Measuring Tape | 67019 | This model uses compression bands, which makes it easier to identify how firmly the tape measure is being pulled. The compression band is not necessary, but highly recommended. |
| MUAC measuring tape (ShorrTape© Measuring Tape) | Weigh and Measure, LLC | STape | The tape measure should be flexible with a single or double slotted insertion window. |