

ECG mode and signal filtering.

| | |
|------------------------------------|---------------------------|
| mode | ECG ; RR-only |
| filter mode | auto, according to set HR |
| expected HR | bpm > 300 |
| baseline removal filter width (ms) | 100.00 |
| noise removal filter width (ms) | 1.00 |
| notch filter | 50 Hz |
| spike removal filter | off |
| drop-out detection mode | off |
| max RR lengths (ms) | 900.00 |
| RR from adjusted R peaks | off |
| RR_only settings mode | Xsmall: mouse |
| R peak width (ms) | 10.00 |
| PR width (ms) | 20.00 |
| RT width (ms) | 50.00 |
| max inter beat artefact (%) | 50.00 |
| R to other amplitude ratio | 3.00 |
| R peak sign | positive |
| compute extra parameter | off |

help

print

cancel

ok