



Run

Run name xmipp3 - preprocess micrograp **Comment**

Run mode Continue Restart **Host** localhost

Parallel Threads MPI **Use queue?** Yes No

Wait for

Expert Level Normal Advanced

Preprocess

Preprocess

Input micrographs xmipp3 - movie alignment CRR+OF.outputMicrographs

Operations are performed in the order shown below

Crop borders? Yes No

Take logarithm? Yes No

Remove bad pixels? Yes No

Multiple of Stddev

Invert contrast? Yes No

Downsample micrographs? Yes No

Downsampling factor

Denoising Yes No

Gaussian filter Yes No

Highpass filter Yes No

Lowpass filter Yes No

Normalize micrograph? Yes No

Close

Save

Execute