

EPU EPU

Preparation Auto Functions Atlas EPU

Start Stop Presets Thon Ring


Execution Preset Selection

Calibrations


Auto-Functions (TEM)

- Autofocus
- Auto-eucentric by beam tilt
- Auto-eucentric by stage tilt
- Autostigmatate**
- Autocoma
- Drift stabilization
- Auto Zero-Loss

Image



Fit Image



Messages

0 Errors 1 Notifications

The exposure time was revised by the camera to 2.36 seconds.

Status

- Starting Autostigmatate...
- CTF estimation started...
- Calculating power spectrum...
- Estimating astigmatism...
- Elliptical averaging...
- Baselining radial profile...
- Estimating defocus...
- CTF estimation completed.
- Phase shift 0.00π
- Defocus $-1.100 \mu\text{m}$
- Astigmatism $3 \text{ nm } (88.2^\circ)$

Histogram

Auto Filter

Image Information

Athena Login