

96-well plate map schematic for ADCC assay

Plate	1	2	3	4	5	6	7	8	9	10	11	12
A	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS
B	PBS	Target only	Target + Dara	Target + NK	Target + NK+ Dara	Target + NK+ Iso	PBS	PBS	PBS	PBS	PBS	PBS
C	PBS	Target only	Target + Dara	Target + NK	Target + NK+ Dara	Target + NK+ Iso	PBS	PBS	PBS	PBS	PBS	PBS
D	PBS	Target only	Target + Dara	Target + NK	Target + NK+ Dara	Target + NK+ Iso	PBS	PBS	PBS	PBS	PBS	PBS
E	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS
F	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS
G	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS
H	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS	PBS
	Target cells	25 K	25 K	25 K	25 K	25 K						
	NK cells	0 K	0 K	12.5 K	12.5 K	12.5 K						
	Antibody	-	Dara	-	Dara	Iso						