

Materials List for:

Preparation of Viral DNA from Nucleocapsids

Moriah L. Szpara¹, Yolanda R. Tafuri¹, L. W. Enquist¹

¹Department of Molecular Biology, Princeton University

Correspondence to: L. W. Enquist at lenquist@Princeton.edu

URL: <https://www.jove.com/video/3151>

DOI: [doi:10.3791/3151](https://doi.org/10.3791/3151)

Materials

Name	Company	Catalog Number	Comments
Freon (1,1,2-trichloro-1,2,2-trifluoroethane)	Fisher Scientific	T178-4	Check with your institution for guidelines on appropriate disposal of Freon-containing waste, or see Mendez et al. for potential Freon alternatives. ⁹
Phase lock gel tubes, Heavy 15ml capacity	5 PRIME	2302850	Optional
Polyallomer ultracentrifuge tubes	Beckman Coulter Inc.	331372*	*Select tubes appropriate for your own ultracentrifuge; these are included as an example only
NP-40 / IGEPAL	Sigma-Aldrich	I-3021	
PK-15 cells	ATCC	CCL-33	
Vero cells	ATCC	CCL-81	
PBS	Hyclone	SH30028.03	