

Materials List for

A Versatile, Behavioral Method to Investigate Thyroid Hormone Effects on Cerebellar Function

Ayane Ninomiya¹, Izuki Amano¹, Noriyuki Koibuchi¹

¹Department of Integrative Physiology, Gunma University Graduate School of Medicine

Corresponding Author

Noriyuki Koibuchi

nkoibuch@gunma-u.ac.jp

Citation

Ninomiya, A., Amano, I., Koibuchi, N. A Versatile, Behavioral Method to Investigate Thyroid Hormone Effects on Cerebellar Function. *J. Vis. Exp.* (200), e65940, doi:10.3791/65940 (2023).

Date Published

October 6, 2023

DOI

10.3791/65940

URL

jove.com/video/65940

Materials

Name	Company	Catalog Number	Comments
Air puff	DAISO		
Aluminum sticks	CAINZ		2 mm diameter, number of sticks may vary depending on the ladder design. Approximately 30 sticks may be required to build the horizontal ladder (4 cm interval).
Blutack	Bostik		
Plexiglass	CAINZ		110 cm x 20 cm, 110 cm x 10 cm, 2 parts each
Screws	CAINZ		