

Materials List for:

Assaying β -amyloid Toxicity using a Transgenic *C. elegans* Model

Vishantie Dostal¹, Christopher D. Link^{1,2}

¹Institute for Behavioral Genetics, University of Colorado

²Integrative Physiology, University of Colorado

Correspondence to: Christopher D. Link at linkc@colorado.edu

URL: <https://www.jove.com/video/2252>

DOI: [doi:10.3791/2252](https://doi.org/10.3791/2252)

Materials

Name	Company	Catalog Number	Comments
NaCl (6g)	Sigma-Aldrich	S9888	
Agar (40g)	Sigma-Aldrich	A7002	
Peptone (5g)	Sigma-Aldrich	P0556	
ddH ₂ O (2L)			
1M CaCl ₂ (2mL)			147.02mg/mL in ddH ₂ O
Uracil (2mL)			2mg/mL in ddH ₂ O
Cholesterol (2mL)			5mg/mL in ethanol (do not autoclave)
1M MgSO ₄ (2mL)			246.5mg/mL in ddH ₂ O
1M KPO ₄ (50mL)			98mg KH ₂ PO ₄ .mL, 48mg K ₂ HPO ₄ /mL in ddH ₂ O, pH
Erlenmeyer Flask	VWR international		