

Materials List for:

A 1.5 Hour Procedure for Identification of *Enterococcus* Species Directly from Blood Cultures

Margie A. Morgan¹, Elizabeth Marlowe², Susan Novak-Weekly², J.M. Miller², T.M. Painter³, Hossein Salimnia³, Benjamin Crystal⁴

¹Department of Pathology and Laboratory Medicine, Cedars-Sinai Medical Center

²Pasadena, CA, Southern California Permanente Medical Group

³Detroit, Detroit Medical Center

⁴Woburn, MA, AdvanDx

Correspondence to: Margie A. Morgan at morganma@cshs.org

URL: <https://www.jove.com/video/2616>

DOI: [doi:10.3791/2616](https://doi.org/10.3791/2616)

Materials

Reagents

E. faecalis- PNA FISH is comprised of the following kit components:

- Fixation Solution Fixation Solution 3 mL phosphate-buffered saline with detergent.
- *E. faecalis*- PNA *E. faecalis*- PNA 1.5 mL PNA probes in hybridization. solution. Contains 30% formamide.
- 60x Wash Solution 60x Wash Solution 50 mL Tris-buffered saline with detergent
- Mounting Medium Mounting Medium 3 mL photobleaching inhibitor in Glycerol
- Quality control slides
- [AdvanDX](#)