

Materials List for:

Induction and Monitoring of Adoptive Delayed-Type Hypersensitivity in Rats

Christine Beeton¹, K. George Chandy¹

¹Department of Physiology and Biophysics, University of California, Irvine (UCI)

Correspondence to: Christine Beeton at cbeeton@uci.edu, K. George Chandy at gchandy@uci.edu

URL: <https://www.jove.com/video/325>

DOI: [doi:10.3791/325](https://doi.org/10.3791/325)

Materials

Name	Company	Catalog Number	Comments
DMEM	GIBCO, by Life Technologies	12430-062	500ml
Sodium pyruvate	Sigma-Aldrich	S8636	1 mM final, in Stimulation medium
Penicillin – Streptomycin – L- Glutamine	Cambrex/BioWhittaker	17-718R	100 u/ml - 100 µg/ml - 4 mM respectively, final, in Stimulation Medium
Non essential amino acids	Sigma-Aldrich	R7131	NEAA, 1% final, in Stimulation Medium
RPMI vitamins	Sigma-Aldrich	R7256	1% final, in Stimulation medium
beta-mercapt–thanol	Sigma-Aldrich	M6250	2ME, 50 µM final, in Stimulation Medium
Syngeneic rat serum			1% final. See explanation in Protocol.
Ovalbumin	Sigma-Aldrich	A5503	10 µg/ml final, in Stimulation Medium.