

Materials List for:

Determining the Phagocytic Activity of Clinical Antibody Samples

Elizabeth G. McAndrew¹, Anne-Sophie Dugast¹, Anna F. Licht¹, Justin R. Eusebio¹, Galit Alter¹, Margaret E. Ackerman²

¹Massachusetts General Hospital, Ragon Institute of MGH, MIT, and Harvard

²Thayer School of Engineering, Dartmouth College

Correspondence to: Margaret E. Ackerman at margaret.e.ackerman@dartmouth.edu

URL: <https://www.jove.com/video/3588>

DOI: [doi:10.3791/3588](https://doi.org/10.3791/3588)

Materials

Name	Company	Catalog Number	Comments
EZ-Link Micro Sulfo-NHS-LC-Biotinylation Kit	Thermo Fisher Scientific, Inc.	21935	
Amicon Ultra-15 Centrifugal Filter Unit with Ultracel-50 membrane	EMD Millipore	UFC905024	Filter unit size may vary according to antigen size
FluoSpheres NeutrAvidin labeled microspheres, 1.0 µm, yellow-green fluorescent (505/515)	Invitrogen	F-8776	Manufacturers produce this product in various colors
Melon Gel IgG Purification Kit	Thermo Fisher Scientific, Inc.	45212	Care should be taken to verify the purity of Ab samples by SDS-page.