

Materials List for:

Use of a Robot for High-throughput Crystallization of Membrane Proteins in Lipidic Mesophases

Dianfan Li¹, Coilin Boland¹, Kilian Walsh¹, Martin Caffrey¹

¹Membrane Structural and Functional Biology Group, Schools of Medicine and Biochemistry & Immunology, Trinity College Dublin

Correspondence to: Martin Caffrey at martin.caffrey@tcd.ie

URL: <https://www.jove.com/video/4000>

DOI: [doi:10.3791/4000](https://doi.org/10.3791/4000)

Materials

Name	Company	Catalog Number	Comments
Brayer (roller)	Fisher Scientific	50820937	Tool
Gas-tight syringes	Hamilton	81030	Tool
Glass coverslips	Marienfeld	01029990911	Disposable
Glass plates	Marienfeld	1527127092	Disposable
Gryphon LCP Robot	Art Robbins	NA	Tool
<i>In meso</i> robot	Anachem/Gilson	NA	Tool
Lab notebook	Various	NA	Tool
Mosquito LCP Robot	TTP Labtech	NA	Tool
Perforated double-stick spacer tape	Saunders Corporation (hole-punched)	customized	Disposable
Precipitant solutions	Various	Various	Reagent
Purified water	Millipore	NA	Reagent
Rain-X	Shell Car Care	80199200	Reagent
Syringe tips	Hamilton	7770-020 (gauge 22)	Tool
Tissues	Various	Various	Disposable
Water bottle	various	NA	Reagent