

Materials List for:

Passaging HuES Human Embryonic Stem Cell-lines with Trypsin.

Erin Trish¹, John Dimos¹, Kevin Eggan¹

¹Department of Molecular and Cell Biology, Harvard

Correspondence to: Erin Trish at etrish@mcb.harvard.edu, John Dimos at jdimos@fas.harvard.edu, Kevin Eggan at eggan@mcb.harvard.edu

URL: <https://www.jove.com/video/49>

DOI: [doi:10.3791/49](https://doi.org/10.3791/49)

Materials

Name	Company	Catalog Number	Comments
HuES media 651.5 ml			
KO-DMEM 500 ml			
PenStrep 6.5 ml			
GlutaMAX TM 6.5 ml			
NEAA 6.5 ml			
2-mercapt-ethanol 0.5ul			
KO Serum Replacement 65 ml			
Plasmanate 65 ml			
bFGF2 (10 ng/mL final) 0.65 ml			