

Materials List for:

A Procedure to Study the Effect of Prolonged Food Restriction on Heroin Seeking in Abstinent Rats

Firas Sedki¹, Tracey D'Cunha¹, Uri Shalev¹

¹Department of Psychology, Center for Studies in Behavioural Neurobiology/Groupe de Recherche en Neurobiologie Comportementale, Concordia University

Correspondence to: Uri Shalev at uri.shalev@concordia.ca

URL: <https://www.jove.com/video/50751>

DOI: [doi:10.3791/50751](https://doi.org/10.3791/50751)

Materials

Name	Company	Catalog Number	Comments
Name	Company	Catalogue Number	
Ketamine (Vetalar, 100 mg/ml)	Bioniche Animal Health	Obtained through local distributor	
Xylazine (Rompun, 100 mg/ml)	Bayer	Obtained through local distributor	
Heroin HCl	National Institute for Drug Abuse, Research Triangle Park, NC, USA		
Anafen injection 1 mg/ml Vial/20 ml	CDMV, Canada	12868	
Equipment			
Name	Company	Catalogue Number	
Silastic tubing (ID 0.02, OD 0.037)	Fisher Scientific	1118915A	
GE marine silicon	GE	SE-1134	
Tygon tubing (ID 0.02, OD 0.060)	VWR	63018-044	
Cannulae (22 G, 5-up)	Plastics One	C313G-5up	
Liquid swivels, plastic, 22 G	Lomir Biomedical, Inc.	RSP1	
Fixed speed infusion pump (3.3 rpm)	Coulbourn Instruments	A73-01-3.3	
Rat test cage	Coulbourn Instruments	H10-11R-TC	
Stainless steel grid floor	Coulbourn Instruments	H10-11R-TC-SF	
House light-rat	Coulbourn Instruments	H11-01R	
Single high-bright cue-rat	Coulbourn Instruments	H11-03R	
Tone module 2.9 kHz	Coulbourn Instruments	H12-02R-2.9	
Retractable lever-rat	Coulbourn Instruments	H23-17RA	
Balance arm	Coulbourn Instruments	H29-01	
Sound attenuation boxes	Concordia University	Home made	
System controller 2	Coulbourn Instruments	SYS CTRL 2	
System power base	Coulbourn Instruments	H01-01	
Habitest Universal Linc	Coulbourn Instruments	H02-08	
Environment connection board & Linc cable	Coulbourn Instruments	H03-04	
Graphic State Notation 3	Coulbourn Instruments	GS3	