

Materials List for:

# Methods for Cell-attached Capacitance Measurements in Mouse Adrenal Chromaffin Cell

Kelly T. Varga<sup>1</sup>, Zhongjiao Jiang<sup>1</sup>, Liang-Wei Gong<sup>1</sup>

<sup>1</sup>Department of Biological Sciences, University of Illinois at Chicago

Correspondence to: Liang-Wei Gong at [lwong@uic.edu](mailto:lwong@uic.edu)

URL: <https://www.jove.com/video/52024>

DOI: [doi:10.3791/52024](https://doi.org/10.3791/52024)

## Materials

Name	Company	Catalog Number	Comments
Poly-D-Lysine	Sigma	P0899	
DMEM		15066024	Keep out of UV
Dulbecco's Modified Eagle Medium	Life Technologies		
Cover Glass	Carolina Biological	633029	12 mm
Penicillin Streptomycin	Life Technologies	15140122	100 ml
Insulin-Trans-Sel-X	Life Technologies	51500056	Only thaw on ICE!
Papain	Worthington	39S11614	
EPC-7 plus patch amplifier	HEKA		
BNC-2090 data acquisition board	National Instruments		
Igor data acquisition software	Wavemetrics		
P-97 pipette puller	Sutter Instruments		
Microforge	Scientific Instruments		
Borosilicate glass capillaries	Sutter Instruments	B150-110-10	Outer diameter: 1.5 mm Inner diameter: 1.10 mm Length: 10 cm