

Materials List for:

Three Different Protocols of Corneal Collagen Crosslinking in Keratoconus: Conventional, Accelerated and Iontophoresis

Nacim Bouheraoua^{*1,2,3,4}, Lea Jouve^{*1,2,3,4}, Vincent Borderie^{1,2,3,4}, Laurent Laroche^{1,2,3,4}

¹Quinze-Vingts National Ophthalmology Hospital

²INSERM UMR S 968, Institut de la Vision

³Sorbonne Universités, UPMC Univ Paris 06

⁴CNRS, UMR 7210

* These authors contributed equally

Correspondence to: Nacim Bouheraoua at nacim.bouheraoua@gmail.com

URL: <https://www.jove.com/video/53119>

DOI: [doi:10.3791/53119](https://doi.org/10.3791/53119)

Materials

Name	Company	Catalog Number	Comments
		Riboflavin Product number	
C-CXL	Sooft SPA, Montegiorgio, Italy	Ricrolin 468465-6	
A-CXL	Avedro Inc, Waltham, Massachusetts	VibeX 520-01863-006	
I-CXL	Sooft SPA, Montegiorgio, Italy	Ricrolin+ 975481-6	Passive electrode: PROTENS ELITE 4848LE/ Active electrode: IONTOFOR CXL
		UVA Machine	
X-Vega			UVA: 3 mW/cm ² 30 min
KXL System			UVA: 30 mW/cm ² 10 min
X-Vega			UVA: 10 mW/cm ² 9 min