

Materials List for:

Repeated Blood Collection for Blood Tests in Adult Zebrafish

Liqing Zang¹, Yasuhito Shimada², Yuhei Nishimura², Toshio Tanaka², Norihiro Nishimura¹

¹Department of Translational Medical Science, Graduate School of Medicine, Mie University

²Department of Molecular and Cellular Pharmacology, Pharmacogenomics, and Pharmacoinformatics, Graduate School of Medicine, Mie University

Correspondence to: Liqing Zang at liqing@doc.medic.mie-u.ac.jp

URL: <https://www.jove.com/video/53272>

DOI: [doi:10.3791/53272](https://doi.org/10.3791/53272)

Materials

Name	Company	Catalog Number	Comments
Glass capillaries with filament	Narishige	GD-1	1.0 mm outer diameter.
Needle puller	Narishige	PC-10	To produce the needles
Heparin	Wako Pure Chemical Industries	081-00136	For heparinization
Aspirator tube assembly	Drummond	2-040-000	For blood collection
Bulb dispenser	Drummond	1-000-9000	For blood collection
2-phenoxyethanol	Wako Pure Chemical Industries	163-12075	For anesthetizing the fish
DRI-CHEM3500V	Fujifilm	-	For hemoglobin measurement
DRI-CHEM Slides	Fujifilm	Hb-WII	For hemoglobin measurement
Glutest Neo Super	Sanwa Kagaku Kenkyusho	-	For blood glucose measurement
Wako L-type TG kit	Wako Pure Chemical Industries	464-44201	For TG measurement
Wako L-type CHO kit	Wako Pure Chemical Industries	460-44301	For total cholesterol measurement
Parafilm M	Alcan Packaging	PM996	To expel the blood on