

Materials List for:

Isolation of Viral Replication Compartment-enriched Sub-nuclear Fractions from Adenovirus-infected Normal Human Cells

Paloma Hidalgo^{1,2}, Ramón A. Gonzalez¹

¹Centro de Investigación en Dinámica Celular, Instituto de Investigación en Ciencias Básicas y Aplicadas, Universidad Autónoma del Estado de Morelos

²Instituto de Biotecnología, Universidad Nacional Autónoma de México

Correspondence to: Ramón A. Gonzalez at rgonzalez@uaem.mx

URL: <https://www.jove.com/video/53296>

DOI: [doi:10.3791/53296](https://doi.org/10.3791/53296)

Materials

Name	Company	Catalog Number	Comments
DMEM	Gibco	12100-046	Warm in 37 °C water bath before use
Fetal Bovine Serum	Gibco	12484-028	
Sucrose, Ultra Pure	Research Organics	0928S	Prepare a 2.55 M stock solution and store at 4 °C
Dounce homogenizer	Kontess Glass Company	884900-0000	
Branson 1800 Ultrasonic Bath	Branson	Z769533 SIGMA	Turn on 15 min before use.
Goat anti-Mouse IgG1 Secondary Antibody, Alexa Fluor 488 conjugate	Life Technologies	A-21121	Use at a 1:2,000 dilution in PBS
Silane-Prep Slides	Sigma	S4651-72EA	Open in a laminar flow cabinet
SuperSignal West Pico Chemiluminescent Substrate	Pierce ThermoScientific	34080	