

Materials List for:

Autologous Microfractured and Purified Adipose Tissue for Arthroscopic Management of Osteochondral Lesions of the Talus

Riccardo D'Ambrosi^{1,2}, Cristian Indino¹, Camilla Maccario¹, Luigi Manzi¹, Federico Giuseppe Usuelli¹

¹IRCCS Istituto Ortopedico Galeazzi - C.A.S.C.O. Piede e Caviglia

²Dipartimento di Scienze Biomediche per la Salute, Università degli Studi di Milano

Correspondence to: Riccardo D'Ambrosi at riccardo.dambrosi@hotmail.it

URL: <https://www.jove.com/video/56395>

DOI: [doi:10.3791/56395](https://doi.org/10.3791/56395)

Materials

Name	Company	Catalog Number	Comments
PROCESS KIT - PROCESSING KIT FOR FAT TISSUE	LIPOGEMS	LG PK 60	Lipogems Kit to obtain microfractured and purified ADSCs
HINTERMANN SPREADER	INTEGRA	119654	The spreader allow to access most of the talar dome, in particular in case of posterior lesion
CUP CURETTE	ARTHREX	AR-8655-02	To remove the damaged cartilage and necrotic and sclerotic bone
CHONDRAL PICK 30° TIP / 60° TIP	ARTHREX	AR-8655-05 AR-8655-06	To perform microperforation at the level of the lesion, stimulating bleeding from the subchondral bone
SHAVER	ARTHREX	AR-7300SR	To clean the joint and aspirate water