

Materials List for:

Recording Horizontal Saccade Performances Accurately in Neurological Patients Using Electro-oculogram

Yasuo Terao¹, Hideki Fukuda², Yusuke Sugiyama³, Satomi Inomata-Terada¹, Shin-ichi Tokushige⁴, Yoshikazu Ugawa⁵

¹Department of Cell Physiology, Kyorin University

²Segawa Memorial Neurological Clinic for Children

³Department of Neurology, University of Tokyo

⁴Department of Neurology, Kyorin University

⁵Department of Neurology, Fukushima Medical University

Correspondence to: Yasuo Terao at yterao-ky@umin.ac.jp

URL: <https://www.jove.com/video/56934>

DOI: [doi:10.3791/56934](https://doi.org/10.3791/56934)

Materials

Name	Company	Catalog Number	Comments
Electrode	Nihon-Kohden (Tokyo, Japan)	NS111-115	cup electrode
Electrode paste	Nihon-Kohden (Tokyo, Japan)	Gelaid Z-101BA	gel electrode paste to fill in the cup electrode
Adhesive tape	Nihon-Kohden (Tokyo, Japan)	H261	double-stick tape for fixating the electrode
DC-amplifier	Nihon-Kohden (Tokyo, Japan)	AN-601G	amplifier for EOG
video-based eye tracking system	SR research (Mississauga, Ontario, Canada)	Eyelink II	eye tracking system for recording VOG
Filter	NF corporation	MS-521	filter for the EOG signal