

**Materials List for:**

**Measurement of Spatial Stability in Precision Grip**

Ryoto Teshima<sup>1</sup>, Naoto Noguchi<sup>2</sup>, Ryota Fujii<sup>3</sup>, Ken Kondo<sup>4</sup>, Koji Tanaka<sup>4</sup>, Bumsuk Lee<sup>4</sup>

<sup>1</sup>Division of Rehabilitation Service, Geriatrics Research Institute and Hospital

<sup>2</sup>Faculty of Rehabilitation, Gunma University of Health and Welfare

<sup>3</sup>Department of Rehabilitation, Gunma Chuo Hospital

<sup>4</sup>Graduate School of Health Sciences, Gunma University

Correspondence to: Bumsuk Lee at [leebumsuk@gunma-u.ac.jp](mailto:leebumsuk@gunma-u.ac.jp)

URL: <https://www.jove.com/video/59699>

DOI: [doi:10.3791/59699](https://doi.org/10.3791/59699)

**Materials**

Name	Company	Catalog Number	Comments
Alcohol swab			Wipe participant's finger pulps
Compressor	Nitta Corporation		Apply pressure to the sensor seats
Computer			
Controller of compressor	Nitta Corporation		Use to manipulate the compressor
Double-sides tapes			Use to attach the sensorseats to the iron cube
Iron cube			150-250g, 30×30×30 mm
Sensor connector			Connect the sensorseats to computer.
Sensor sheet			Pressure Mapping Sensor 5027, Tekscan, South Boston, MA, 50 USA
Setting stand			Set the iron cube on it during the measurement
Software; I-SCAN 5027, Ver. 7.51	Nitta Corporation		
Table			Use for the measurement