

Materials List for:

Tuning Degradation to Achieve Specific and Efficient Protein Depletion

J. David Barrass¹, Gonzalo I. Mendoza-Ochoa^{1,2}, Isabella E. Maudlin^{1,3}, Emanuela Sani¹, Jean D. Beggs¹

¹Wellcome Centre for Cell Biology, School of Biological Sciences, University of Edinburgh

²Department of Plant Sciences, University of Cambridge

³Sir William Dunn School of Pathology, University of Oxford

Correspondence to: J. David Barrass at david.barrass@ed.ac.uk

URL: <https://www.jove.com/video/59874>

DOI: [doi:10.3791/59874](https://doi.org/10.3791/59874)

Materials

Name	Company	Catalog Number	Comments
Adenine sulphate	Formedium	DOC0230	
Agar	Formedium	AGA03	
B-estradiol	Sigma Aldrich	E2758-1G	10 mM solution in ethanol. Store at -20 °C
DMSO	Alfa Aesar	42780	DMSO should be solid at 4 °C
Glucose	Fisher Scientific	G/0500/60	
IAA 1H-Indole-3-acetic acid	Across Organics	122150100	Auxin analogue. 1.5 M in DMSO. The solution will be a russet colour and darken as time goes on; a deep red solution should be discarded and a new one made. Store at -20 °C.
Methanol	Fisher Scientific	M/4000/PC17	CAUTION Toxic and flammable
Phusion High-Fidelity DNA Polymerase	NEB	M0530	
Peptone	Formedium	PEP03	
SCSM single drop-out –ura	Formedium	DSCS101	
Yeast Extract	Formedium	YEA03	
Yeast nitrogen base without amino acids with amonium sulphate	Formedium	CYN0410	