

Materials List for:

A Method for Tracking the Time Evolution of Steady-State Evoked Potentials

Pavel Prado-Gutiérrez¹, Mónica Otero^{1,2}, Eduardo Martínez-Montes³, Alejandro Weinstein^{1,4}, María-José Escobar^{1,2}, Wael El-Deredy^{1,4}, Matías Zañartu^{1,2}

¹Advanced Center for Electrical and Electronic Engineering (AC3E), Universidad Técnica Federico Santa María

²Department of Electronic Engineering, Universidad Técnica Federico Santa María

³Neuroinformatics Department, Cuban Neuroscience Center

⁴Centro de Investigación y Desarrollo en Ingeniería, Universidad de Valparaíso

Correspondence to: Pavel Prado-Gutiérrez at pavel.prado@gmail.com

URL: <https://www.jove.com/video/59898>

DOI: [doi:10.3791/59898](https://doi.org/10.3791/59898)

Materials

Name	Company	Catalog Number	Comments
Active electrodes	Biosemi	P32-1020-32ACMS (ABC)	for channels 1-32
Active electrodes	Biosemi	P32-1020-32A (ABC)	for channels 33-64
Active electrodes	Biosemi	8 x TP FLAT	external electrodes
Active-Two acquisition system	Biosemi	version 7.0	EEG acquisition system
alcohol	Salcobrand	Code: 3309011	for cleaning the scalp
Electrode cap 64 channels	Biosemi	CAP MS xx yy	cap
Electrode cap 64 channels	Biosemi	CAP ML xx yy	cap
gel	Biosemi	SIGNA BOX12	conductive gel
Laptop	Asus	Core i7 1TB DD + 128GB SSD 8GB RAM	computer for stimulation
Laptop	Asus	Core i7 1TB DD + 128GB SSD 8GB RAM	computer for recording
LED screen	in-house production	-	The screen consists of four light-emitting diodes (LEDs) situated on the center of a 50x50 cm black screen, as vertexes of a square of 5x5 cm
sterile gauze	Salcobrand	Code: 8730277	for cleaning the scalp