

Materials List for:

Generation of Stable Transgenic *C. elegans* Using Microinjection

Laura A. Berkowitz¹, Adam L. Knight¹, Guy A. Caldwell¹, Kim A. Caldwell¹

¹Department of Biological Sciences, University of Alabama

URL: <https://www.jove.com/video/833>

DOI: [doi:10.3791/833](https://doi.org/10.3791/833)

Materials

Name	Type	Company	Catalog Number	Comments
Agarose	Ultrapure	Invitrogen	15510-027	
Halocarbon Oil, Voltalef	Hulle 10S	elfatochem, France		
Coverglass 18x18 mm		Fisher Scientific	12-548-A	
Coverglass 20x30 mm		Fisher Scientific	12-548-5A	
Coverglass 22x50 mm		Fisher Scientific	12-545-E	
100 ul Capillary		VWR international	53432-921	
Glass Capillary for Needles	"Kwik-Fil"	World Precision Instruments, Inc.	1B100F-4	
Needle Puller	Tool	Narishige International	Model PP-830	
microINJECTOR System		Tritech Research, Inc.	MINJ-1000	Scope, Stage, Manipulator
Dissecting, with Fluorescence	Microscope	Nikon Instruments	SMZ800	
Dissecting	Microscope	Nikon Instruments	SMZ645	