

Materials List for:

Selective Depletion of Microglia from Cerebellar Granule Cell Cultures Using L-leucine Methyl Ester

Joseph Jebelli^{*1}, Thomas Piers^{*2}, Jennifer Pocock³

¹Department of Neurology, University of Washington

²Therapeutic Innovation Group, School of Life and Medical Sciences, University College London

³Department of Neuroinflammation, University College London

*These authors contributed equally

Correspondence to: Joseph Jebelli at jjebelli@uw.edu

URL: <https://www.jove.com/video/52983>

DOI: [doi:10.3791/52983](https://doi.org/10.3791/52983)

Materials

Name	Company	Catalog Number	Comments
Forceps	Sigma-Aldrich	F4142	The curved end facilitates removal of the cerebellum
Micro-dissecting scissors	Sigma-Aldrich	S3146	Straight, sharp point facilitates rodent P4-7 dissection
L-leucine methyl ester hydrochloride	Sigma-Aldrich	7517-19-3	
EBSS solution	Sigma-Aldrich	E7510-500 ml	
Poly-D-lysine	Sigma-Aldrich	27964-99-4	Coat coverslips 1 day before use
Bovine serum albumin (BSA)	Sigma-Aldrich	A9418	
Phosphate buffered saline (PBS)	Sigma-Aldrich	P4417	
DNase	Sigma-Aldrich	D5025	
Soybean trypsin inhibitor	Sigma-Aldrich	T6414	
Mouse anti-ED1 antibody	Abcam	ab31630	