

Materials List for:

In Vitro Assay to Study Tumor-Macrophage Interaction

Zhenyi An¹, William A. Weiss^{1,2,3}

¹Department of Neurology, University of California, San Francisco

²Helen Diller Family Comprehensive Cancer Center, University of California, San Francisco

³Departments of Pediatrics and Neurological Surgery, University of California, San Francisco

Correspondence to: Zhenyi An at anzhenyi@gmail.com

URL: <https://www.jove.com/video/59907>

DOI: [doi:10.3791/59907](https://doi.org/10.3791/59907)

Materials

Name	Company	Catalog Number	Comments
0.1 µm filtration cup	Thermo fisher	566-0010	
0.45 µm filter unit	Millipore	SLHA033SS	
10 mL serological pipettes	Olympus plastics	12-104	
15mL sterile centrifuge tubes	Olympus plastics	28-103	
1 mL pipette tip	ART molecular bioproducts	2779-RI	
2 mL aspirating pipet	Falcon	357558	
24-well plate	Millipore	ECM507	Part of ECM507, or can be purchased separately
4x lysis buffer	Millipore	ECM507	Part of ECM507, or can be purchased separately
5 µm Transwell insert	Millipore	ECM507	Part of ECM507, or can be purchased separately
75cm ² flask	Corning	430641U	
Accutase	Innovative cell technologies	AT-104	
B27	Gibco	12587-010	
CyQuant Dye	Millipore	ECM507	Part of ECM507, or can be purchased separately
DMEM	Gibco	11965-092	
DMEM:F12	Gibco	10565-018	
EGF	Peptrotech	AF-100-15	
FBS	Gibco	26140	
FGF	Peptrotech	100-18B	
IMDM	Gibco	12440-053	
PBS	Gibco	14190-144	
Pen Strep	Gibco	15140-122	
Trypan blue	Biorad	1450021	