## Materials List for: Training Dogs for Awake, Unrestrained Functional Magnetic Resonance Imaging

Lily R. Strassberg\*<sup>1</sup>, L. Paul Waggoner\*<sup>2</sup>, Gopikrishna Deshpande\*<sup>1,3,4,5</sup>, Jeffrey S. Katz\*<sup>1,3,4,5</sup>

<sup>1</sup>Department of Psychology, Auburn University

<sup>2</sup>Canine Performance Sciences Program, College of Veterinary Medicine, Auburn University

<sup>3</sup>AU MRI Research Center, Department of Electrical & Computer Engineering, Auburn University

<sup>4</sup>Alabama Advanced Imaging Consortium

<sup>5</sup>Center for Neuroscience, Auburn University

\*These authors contributed equally

Correspondence to: Jeffrey S. Katz at katzjef@auburn.edu

URL: https://www.jove.com/video/60192 DOI: doi:10.3791/60192

## **Materials**

Name	Company	Catalog Number	Comments
Acrylic Mock Radiofrequency Coil	Menards	TU59018594	Mock Radiofrequency (RF) Coil: 8" diameter x 4' Concrete Form Tube. Makes four mock RF coils; cut form tube in four even lengths for four 8" diameter x 1' mock RF coils.
Agility Tunnel	J&J Dog Supplies	TT053	Open Agility Training Tunnel
Bluetooth Speaker	Sharkk	SP-SK896WTR-GRY	Portable Scan Audio Playback: Waterproof Bluetooth Speaker Sharkk 20 IP67 Bluetooth Speaker Outdoor Pool Beach and Shower Portable Wireless Speaker
Cardboard Concrete Form Tube	Menards	TU10120014	Stationary Mock MRI Bore: Sonotube 24" diameter x 12' Standard Wall Water-Resistant Concrete Form. Makes two mock bores; cut form tube in half for two 24" diameter x 6' bores.
Chuckit Ball	Chuckit!	17030	Toy Reward: Chuckit! Ultra Ball
Decibel X	Skypaw		Decibel meter phone app
Exercise Mat			Foam chin rest: cut mat in half lengthwise. Roll up, and secure roll with hot glue. Cut chin-size notch in center with X-ACTO knife. Hot- glue velcro to bottom surface.
Folding Table			3' x 6' folding table
Microfiber Car Wax Applicator Pad	Viking Car Care	862400	Viking Car Care Microfiber Applicator Pads
Natural Balance Treat Log	Natual Balance	236020	Food Reward: E.g., Chicken Formula Dog Food Roll, 3.5-lb roll
Plywood			Platform: 2"x4"x6' length of wood affixed to 3'x6' plywood board. Hot glue exercise mat on plywood board for traction. Braces: 3 4x4x4" cubes cut at 45-degree angle affixed to ends of 1"x4"x3' lengths of wood. Makes 3 braces.
Sand Bags	J&J Dog Supplies	AG155	J&J Professional Quality Sandbags x 2



Speaker System	Pioneer Electrics	Stationary Scan Audio Playback: Pioneer HTD645DV 5 Disk DVD Home Theater System with Wireless Surround Speakers. Operating Instructions.
Towel		standard towel