

Materials List for:

Protocol for RNAi Assays in Adult Mosquitoes (*A. gambiae*)

Lindsey Garver¹, George Dimopoulos¹

¹Malaria Research Institute, Bloomberg School of Public Health, Johns Hopkins University

Correspondence to: George Dimopoulos at gdimopou@jhsph.edu

URL: <https://www.jove.com/video/230>

DOI: [doi:10.3791/230](https://doi.org/10.3791/230)

Materials

Name	Type	Company	Catalog Number	Comments
Megascript	kit	Ambion		for dsRNA synthesis; T7 RNA polymerase-mediated in vitro transcription reaction
RNeasy		Qiagen		For dsRNA purification
dsRNA				should be quantified and adjusted to 3 µg/ µL in water
fine-tipped forceps				
glass slide				
glass microcapillary needles				
microinjector		Drummond Scientific	Nanoject II	
light microscope				mounted above a cold block
Petri dish				covered, glass
paper towels				
10% sucrose	food			
<i>A. gambiae</i>	Animal			mosquitos